

# ACETABULAR REAMER BRUSH

Plastic handle, rigid nylon bristles + double ended

## PRODUCT INFORMATION:

- Designed exclusively for cleaning acetabular bone reamers.
- 7.75" overall length
- Plastic handle, rigid nylon bristles + double ended
- Anti-bacterial medical grade bristles
- The "Pig Tail" rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes.



## ORDERING INFORMATION:

**Description:** Acetabular Reamer Brush

**Product Number:** SH-SUB-08-1

**Packaging:** 3/PK