


# MAGNETIC INSTRUMENT MATS

Non-sterile, single use, not made with natural rubber latex + autoclavable


## PRODUCT INFORMATION:

- Designed to provide a safe way to transfer sharp surgical instruments with the sterile field during procedures.
- Available in 2 sizes  
Small - 16" x 10"  
Large - 16" x 20"
- Strong magnets prevent instruments from falling off the sterile field.

SMALL - SH-MIM-1



LARGE - SH-MIM-2



## ORDERING INFORMATION:

**Description:** Small Magnetic Instrument Mat

**Product Number:** SH-MIM-1

**Packaging:** 24/BX

**Description:** Large Magnetic Instrument Mat

**Product Number:** SH-MIM-2

**Packaging:** 24/BX